



U2™ Knee

All Poly Tibial Component

- Durable, lowered cost, elimination of backside wear, and the ease of removal if necessary. [1,2]
- Multiple articular selections: CR, PS, UC
- Deeper cement recess with rougher surface finishing to enhance cement fixation and anti-rotation
- Featured instrument design ensures safer implantation
- Material: UHMWPE

Reference

[1] All-Polyethylene Versus Metal-Backed Tibial Components—An Analysis of 27,733 Cruciate-Retaining Total Knee Replacements from the Swedish Knee Arthroplasty Register. Asgeir Gudnason, Nils P. Hailer, Annette W-Dahl, Martin Sundberg, Otto Robertsson., *J Bone Joint Surg Am.* 2014

[2] The role of the cemented all-polyethylene tibial component in total knee replacement: a 30-year patient follow-up and review of the literature. Thomas J. Blumenfeld, Richard D. Scott., *Knee.* 2010



Order Information



Thickness Sizes	Cruciate Retaining				Posterior Stabilized				Ultra-Congruent			
	9mm	11mm	13mm	15mm	9mm	11mm	13mm	15mm	9mm	11mm	13mm	15mm
#1	22031011	22031012	22031013	22031014	22031211	22031212	22031213	22031214	22031411	22031412	22031413	22031414
#2	22031021	22031022	22031023	22031024	22031221	22031222	22031223	22031224	22031421	22031422	22031423	22031424
#3	22031031	22031032	22031033	22031034	22031231	22031232	22031233	22031234	22031431	22031432	22031433	22031434
#4	22031041	22031042	22031043	22031044	22031241	22031242	22031243	22031244	22031441	22031442	22031443	22031444
#5	22031051	22031052	22031053	22031054	22031251	22031252	22031253	22031254	22031451	22031452	22031453	22031454
#6	22031061	22031062	22031063	22031064	22031261	22031262	22031263	22031264	22031461	22031462	22031463	22031464
#7	22031071	22031072	22031073	22031074	22031271	22031272	22031273	22031274	22031471	22031472	22031473	22031474



Catalog No. 9304-1103
Description Modular Handle



Catalog No. 9404-5103
Description APT Component Impactor



Each Step We Care



www.uoc.com.tw

Contact Us

■ Taiwan

United Orthopedic Corporation
12F, No. 80, Sec. 1, Chenggong Rd, Yonghe Dist.,
New Taipei City 23452, Taiwan
Tel: +886 2 2929 4567 Fax: +886 2 2922 4567

■ China

United Medical Instrument Co. Ltd.
Room 402, Building 7, No.697, Ling ShiRoad,
Shanghai, China, 200072
Tel: +86 21 3661 8055 Fax: +86 21 6628 3031

■ Japan

United Biomech Japan Inc.
Ginyo Bldg 4th floor, 2-9-40 Kitasaiwai, Nishi-Ku
Yokohama City, Kanagawa, 220-0004, Japan
Tel: +81 45 620 0741 Fax: +81 45 620 0742

■ Europe

UOC (Suisse) SA
Avenue Général Guisan 60A, 1009 Pully, Switzerland.
Tel: +41 (0) 21 634 70 70 Fax: +41 (0) 21 634 70 76

UOC (France)

21, Rue de la Ravinelle, 54000 Nancy, France.
Tel: +33 (0) 3 83 23 39 72 Fax: +33 (0) 3 83 23 39 10

■ EU Representative

mdi Europa GmbH
Langenhagener Strasse 71,
30855 Langenhagen, Germany
Tel: +49 511 3908 9530 Fax: +49 511 3908 9539

■ USA

UOC USA Inc.
20 Fairbanks, Suite 173, Irvine, CA 92618
Tel: +1 949-328-3366 Fax: +1 949-328-3367