

U2[™] MATRIX
Porous Stem

Aim For Perfect Canal Fitting and Filling...

Maximized Fit and Fill

13 different stem configurations provided better fits for a wide varieties of canals



Champagne-flute



Normal



Stovepipe

The U2 Matrix Design


Stem Size (No.)	#7							#7
#6						#6	#6	
#5					#5	#5		
#4				#4	#4			
#3			#3	#3	#3			
#2			#2	#2				
#1		#1						
	Φ9	Φ10	Φ11	Φ12	Φ13	Φ14	Φ15	

Distal Diameter (mm)




New Instrument Available

Easy-to use instrument maximize efficiency in practice



Implant Features:



Teardrop Tip

The teardrop distal end allows easy insertion and minimizes cortical impingement

Stem Configurations

Distal Diameter (mm)	Stem Length (mm)	Proximal Options
Φ9	126	#1
Φ10	132	#2, #3
Φ11	138	#2, #3, #4
Φ12	145	#3, #4, #5
Φ13	152	#5, #6
Φ14	158	#6
Φ15	166	#7

Order Information

Catalog No.	Description (mm)
1104-3212	#1, Φ9
1104-3223	#2, Φ10
1104-3224	#2, Φ11
1104-3233	#3, Φ10
1104-3234	#3, Φ11
1104-3235	#3, Φ12
1104-3244	#4, Φ11
1104-3245	#4, Φ12
1104-3255	#5, Φ12
1104-3256	#5, Φ13
1104-3266	#6, Φ13
1104-3267	#6, Φ14
1104-3278	#7, Φ15



Office

■ Taiwan

United Orthopedic Corporation

12F, No. 80, Sec. 1, Chenggong Rd, Yonghe Dist.,
New Taipei City 23452, Taiwan

Tel: +886 2 2929-4567 Fax: +886 2 2922-4567

■ USA

UOC USA Inc.

20 Fairbanks, Suite 173, Irvine, CA 92618

Tel: +1 949-328-3366 Fax: +1 949-328-3367

■ China

United Medical Instrument Co. Ltd.

Room 3-1-201, No.45, Jiangchang First Road, Shanghai, China

Tel: +86 21 3661-8055 Fax: +86 21 6628-3031

■ EU Representative

mdi Europa GmbH

Langenhagener Strasse 71, 30855 Langenhagen, Germany

Tel: +49-511-3908 9530 Fax: +49-511-3908 9539



www.uoc.com.tw

All Rights Reserved, ©2016 United Orthopedic Corp.

01MP-B-02-03